

TD GLUCO

 **Bluetooth®**

BLOOD GLUCOSE MONITORING SYSTEM



TD-4277B

STRIP
EJECT



No-coding



Bluetooth
pairing



LARGE BACKLIT
SCREEN

LARGE
memory



0.5 µL



5 sec.



4 Alarm settings



7 14 21 28 60 90
Daily average

**GDH-
FAD**

Enzyme technology

KETONES
warning



AST

Alternative
Site
Testing



Accuracy
Safety
Ergonomics



www.tuv.com
ID 0000057751

Measurement adjustable before and after meal



Trust is good, but checking is better!

Features:

- ✓ Bluetooth®-pairing
- ✓ Large backlit lcd- screen
- ✓ Strip- ejection function to prevent blood contact
- ✓ Measurement indication before or after meal (AC/ PC)
- ✓ No- coding
- ✓ Only small blood sample of 0,5 µl required
- ✓ Fast results within 5 seconds

Sample size: 0,5 µL

Response time: 5 seconds

Measuring range: 1,1 to 33,3 mmol/L

Hematocrit range: 35% to 50%



To be used with
TD-4360 test strips

If you have a smartphone, you can use the GlucoCheck app™ to manage the results of your TD-Gluco Bluetooth® blood glucose monitoring system. The link between your meter and your smartphone is via Bluetooth®. The GlucoCheck app™ is suitable for iOS and Android. You can download the correct version via the links below.





BLOOD GLUCOSE MONITORING SYSTEM

Article code	GD100274
Description	TD-GLUCO Bluetooth® starter kit (1) Including 10 test strips
Packaging	Carton
Content	<ul style="list-style-type: none">• TD-GLUCO Bluetooth® meter TD-4277B• TD-GLUCO teststrips (10)• TD-4360 lancing device• Batteries (two EA, 1,5V AAA)• User Manual• Lancets (10)• Protective bag
Carton dimensions (l x w x h)	115 x 55 x 163 mm
Weight	270 grams
EAN code	8718692270817

Article code	GD100222
Description	TD-GLUCO test strips (50) TD-4360
Packaging	Carton
Content	50 test strips manual
Carton dimensions (l x b x h)	44 x 40 x 58 mm
Weight	36 grams
EAN code	8718692270572

The TD-GLUCO Bluetooth® Blood Glucose Monitoring System (TD-4277B) and TD-GLUCO Test Strips (TD-4360) are CE certified medical devices.



Distributor: **GD Medical**
Hastelweg 224, 5652 CL Eindhoven - NL
www.ht-one.nl
Phone number +31(0)40-3031090
info@gdmedical.nl



The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by HT Medical BV is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.